



BANKSIA PRODUCT LIST 2008

FORENSIC DETECTION KITS

<u>Prod. No.</u>	<u>Description</u>	<u>Details</u>
GUN SHOT POWDER RESIDUE KIT		
1800	GSPR Kit (25 tests)	GSPR 1, GSPR 2, GSPR 3 and mittens. Packed in R1 Carrying Case.
Refills		
1810	GSPR 1 (100 cc)	Gun powder residues on hands & near objects
1820	GSPR 2 (100 cc)	Gun powder residues on hands & near objects
1830	GSPR 3 (100 cc)	Gun powder residues on hands & near objects
80300	Mittens for GSPR (25)	Refill of GSPR kit
FINGERPRINT DETECTION		
Large Size Cans		
2067	Nin Plus Ultra (320 cc)	Light colored pourous surfaces
2077	SPR Black (320 cc)	Non-porous, wet light coloured surfaces
2087	SPR White (320 cc)	Non-porous, wet, dark coloured surfaces
2056	Crystal Violet Regular (200 cc)	Adhesive tapes and visualisation of fingerprints treated with cyanoacrilate
Regular Size Cans		
2010	All Print (100 cc)	Light, non-poros, wet or dry (even oily) surfaces
2060	Nin Plus Ultra (100 cc)	Light colored surfaces
2070	SPR Black (100 cc)	Non-porous, wet light colored surfaces
2080	SPR White (100 cc)	Non-porous, wet, dark colored surfaces
Non Aerosol Containers		
2069	Liquidrin (2.5 Litre)	Fingerprints over porous light coloured surfaces
2400	Inkleaner (1 Litre)	Cleaning formula for fingerprint accessories
BIOLOGICAL FLUID RAPID DETECTION		
80FK11	Blood Stain (5 tests)	Multi-device for swab
80FK21	Semen (5 tests)	Device for swab
FORENSIC & MARKING PRODUCTS		
Invisible Marking		
2100	UV Trap (100 cc)	Invisible Objects Spray
2130	UV Crayon	Invisible Objects Marking
2140	UV Lamp	Vision of materials glowing under UV radiation
Transparency spray for envelopes		
1597	See Through (300 cc)	36 Envelopes
1596	See Through (200 cc)	24 Envelopes
1590	See Through (100 cc)	12 Envelopes
Others		
2090	Ferrotrace (100 cc)	Iron Traces Image Detection
2380	Blue (100 cc)	Shoe-Print Developer

Corporate Head Office

3/160 Lytton Road
Bulimba QLD 4171

AUSTRALIA

Ph: +61 7 3902 3000

Fx: +61 7 3217 9869

Free Call: 1300 769944

Email :sales@banksiascientific.com.au

www.banksiascientific.com.au



New Zealand Office

12 Brigade Street
Airport Oaks
Auckland

NEW ZEALAND

Ph: +64 9 921 2626

Fx: +64 9 921 2627

Email: sales@banksiascientific.co.nz

www.banksiascientific.co.nz